Minicall 12

01295.11

Mobiles Rufsystem

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Bild 2: Lautstärkereglер

Please read carefully!

These operating instructions form an integral part of this product and contain important information on how to commission, install and operate it. They must be taken into consideration and given to the end user or other third party when the installer hands over the product. These operating instructions must be retained for future reference.

Damage resulting from abuse or misuse is not covered by warranty. The manufacturer accepts no warranty for consequential losses that may result directly or indirectly from use of this equipment. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The product meets the requirements of the current European and national guidelines for electromagnetic compatibility. Conformity has been established and the relevant statements and documents have been deposited at the manufacturer.

 **It should be noted that for applications that could put human life at risk, by using radio connection there is always present a potential hazard due to extraneous disturbances. The Operating Instructions need to be followed at all times.**

1 Mode of Operation

Minicall12 verifies coded radio frequency signals sent by corresponding transmitters and operates display and signalisation on receipt of a correct signal. The Minicall-12 is not a medical device.

It must be noted that for applications or installations involving wireless transmission and which could have an influence to safety of man or machine there may be a risk of unforeseen operation from extraneous uncontrolled radio signals.

2 Description

Minicall 12 is a self observing calling system with up to 12 calling lines. For each calling line (channel) a maximum of 5 transmitters may be programmed. So, up to 60 transmitters (MT-12 and WT-12) may be used within this system. Once a call is initiated at transmitters side, the Minicall12 will signalize optically and acoustically the corresponding channel and will transmit this call also to portable receivers (R-12). The display of portable receivers will show the channel Number (1..12). A maximum of 5 portable receivers may be programmed.

For the additional display of calls an optional signaler (SG-12) can be used on each channel. In addition to the function as a room signal light, digital outputs can also be used to control external signals.

An active call may be cancelled by pressing the corresponding push button at the Minicall12 or by confirming the call with a portable receiver. If more than 1 call is active, the portable receivers will show them one after the other and only the call appearing in the display may be confirmed by pressing the push button "OK" twice.

If there is to realize a bigger range of operation, up to 5 repeaters may be used.

For any incoming call there is a potential free contact inside the housing which will switch for 1 second (up to 15 seconds selectable). This is for external signalisation devices if needed.

An interruption of radio communication to a transmitter is shown by flashing periodically twice the channel LED and a short BEEP. Just press the push button to stop signalisation if the problem is known.

Low battery of any transmitter is signalized by flashing 5 times periodically the corresponding channel LED.

Summary of signalisation in normal operation mode:

„LED ON" is permanently on	AC-adaptor (Supply) is connected
„LED ON" flashes once a second	Minicall12 is running on battery supply (AC-adaptor missing)
„LED DATA" flashes shortly	signalizes radio communication
„LED DATA" flashes shortly twice	Batteries of Minicall12 are empty or taken off
„Channel LED" (1...12) is permanently on	call (1...12) has been received
„Channel LED" (1...12) flashes periodically twice	radio communication to the corresponding transmitter fails
„Channel LED" (1...12) flashes 5 times periodically	battery capacity of the corresponding transmitter is low

4 Installation and functional check

4.1 Minicall12

Just connect the AC-adaptor and you will hear a starting melody and the green "LED ON" will turn on. The yellow "LED DATA" flashes about twice a second. This will show, that the battery supply is still not activated. Keep this state until later.

4.2 TeachIn (if system is ordered including CG-Service, continue at 4.4)

Transmitter MT-12 / WT-12	Push button
Channel 1	1
Channel 2	2
Channel 3	3
Channel 4	4
Channel 5	5
Channel 6	6
Channel 7	7
Channel 8	8
Channel 9	9
Channel 10	10
Channel 11	11
Channel 12	12

Portable Receiver R-12	Push button
Channel 1	1 + 2
Channel 2	2 + 3
Channel 3	3 + 4
Channel 4	4 + 5
Channel 5	5 + 6

Repeater RT-12	Push button
Channel 1	7 + 8
Channel 2	8 + 9
Channel 3	9 + 10
Channel 4	10 + 11
Channel 5	11 + 12

Signaler SG-12	Push button
Channel 1	1
Channel 2	2
Channel 3	3
Channel 4	4
Channel 5	5
Channel 6	6
Channel 7	7
Channel 8	8
Channel 9	9
Channel 10	10
Channel 11	11
Channel 12	12

Table1: Channel/device togetherness with corresponding push buttons

4.2.1 Teach in transmitter

Press the desired push button for at least 3 seconds until „teach in mode“ is signalized acoustically and the channel LED flashes once a second. Now press the push button „CALL“ of the transmitter that shall control this channel. Therefore keep a distance of about 1 meter. As soon as the transmitters address is learned, the „teach in mode“ will be left and the channel LED will turn off.

In order to turn off the red LED of the transmitter, press the push button „CALL“ once again and then press the channel push button of the Minicall12.

For further transmitters repeat this procedure until all transmitters are taught in. Transmitters that are already known to the Minicall12 will be ignored if you try to reach them in once again. For each channel you may teach in up to 5 transmitters.

If the „teach in mode“ shall be left without any transmitter action, just press once again the corresponding channel push button.

If you pressed the wrong channel LED just press the right one for 3 seconds and the new channel will turn to teach in mode.

If the Minicall12 is set to the „teach in mode“ and no new transmitter is active this mode automatically will finish after a time period of 20 seconds.

4.2.2 Teach in portable receiver R-12

For each receiver channel only one receiver may be used. The teach in process is the same then for transmitters but you have to press two push buttons at the same time (see table 1). For indication of the „teach in mode“ for receivers both LED are flashing (i.e. 1 and 2). The portable receivers automatically send their radio address all 5 seconds. So, if not yet done, switch on the receiver and wait until Minicall12 has learned the desired receiver. The channel used is irrelevant

4.2.3 Teach in repeater RT-12

For each transceiver channel only one transceiver may be used. The teach in process is the same then for transmitters but you have to press two push buttons at the same time (see table 1). For indication of the „teach in mode“ for transceivers both LED are flashing (i.e. 7 and 8).

Switch on the repeater after starting the „teach in mode“. Because of interval duration of automatic repeater activity it could happen, that the „teach in mode“ is left before the repeater is learned. The channel used is irrelevant

4.2.4 Teach in signaler SG-12

For each signaler channel only one signaler may be used. The teach in process is the same then for transmitters (see table 1).

Switch on the signaler after starting the „teach in mode“. Because of interval duration of automatic signaler activity it could happen, that the „teach in mode“ is left before the signaler is learned.

4.3 Erase

4.3.1 Erase all transmitters for one channel

Press the desired push button (i.e. 1) for at least 3 seconds. So the „teach in mode“ is activated. Release the push button and press it again for another 3 seconds until the corresponding LED flashes 5 times and an acoustic signal will confirm deletion of all learned transmitters for this channel. Only all components of the corresponding channel can be deleted.

4.3.2 Erase a portable receiver R-12

Proceed exactly as for erasing a transmitter, but use both push buttons required (i.e. 1 + 2 / see table 1)

4.3.3 Erase a Repeater RT-12

Proceed exactly as for erasing a transmitter, but use both push buttons required (i.e. 7 + 8 / see table 1)

4.3.3 Erase a Signaler SG-12

Proceed exactly as for erasing a transmitter. (see table 1)

4.4 System Check

Distribute all transmitters, signalers (if so) and repeaters (if so) to the place you would like to operate the devices.

Now activate a call at all transmitters by pressing the push button to ensure the correct functioning. The transmitters will turn on their red LED whereas all calls shall be displayed at the Minicall12 by the corresponding red LEDs and the acoustic calling signal sounds all 5 seconds. Additionally the calls must appear in the portable receivers display and signalers light up.

Now press the channel reset push buttons at the Minicall12. Only if all channels are confirmed the calling melody will stop.

Some seconds later the red LED at the transmitters side will turn off and with a short delay also the calling information at the portable receivers will be erased and the lights of the signaler turn off.

If a call of a transmitter placed to far away do not activate the calling channel at the Minicall12, so just place a repeater between this transmitter and the Minicall12. In some cases you may need one more repeater in row to overcome the distance.

Now try again to transmit a call from the transmitter to the Minicall12 in order to confirm the repeaters function.

4.5 Settings (optionally)

4.5.1 Settings Menu Minicall12

You will reach the menu by pressing push button 5 and 6 simultaneously while Power On. This means battery supply and AC-adaptor supply is off and in the moment you switch one of both supplies on push button 5 and 6 have to be pressed.

Minicall12 will show the active settings menu by flashing LED „DATA“ and „ON“ once a second.

With push buttons 1 to 6 you may chose the parameter to be changed and with push buttons 7 to 12 you select the corresponding value (see also following table 2).

Parameter 1: Select one of 6 different acoustic calling signals.

Parameter 2: By default each pressure of a push button is accompanied by a short beep. In case of interruption of the radio communication to the transmitter there will be a beep also. With push button 7 you may switch of this signalization. With push button 9 to12 you may change the repetition frequency of the acoustic calling melody or switch off completely.

Parameter 3: Here the cyclic radio transmission check from the transmitters may be changed. An incoming call is always accepted immediately, but the transmitted reset to the transmitter depends on this chosen time.

Parameter 4: Here the cyclic radio transmission to portable receivers may be changed. Please note that the chosen time may also

pass before a call is transmitted to the receivers.


- Parameter 5: An incoming call closes the internal relay contact as long as this is defined here
- Parameter 6: Do not change, because this is RF frequency select and have to be the same as all parts of the system
- Parameter 7: Here the cyclic radio transmission check from the signalers may be changed. An incoming call is always accepted immediately, but the transmitted reset to the signaler depends on this chosen time.
- Parameter 8: Here the behaviour of the signaler lights (and outputs) may be changed.
- Parameter 9: Here the behaviour of the signaler sound may be changed. This parameter determines whether the sound is activated and at which intervals the tone sounds.
- Parameter 8: You will reach the configuration menu by pressing this push button combination simultaneously while Power On.

Nr.	Push button	Parameter	Value					
			Push button 7	Push button 8	Push button 9	Push button 10	Push button 11	Push button 12
1	1	Calling Melody	Sound 1	Sound 2	Sound 3	Sound 4	Sound 5	Sound 6
2	2	Beep / calling sound repetition	Beep off	Beep on	2 Sec.	5 Sec.	10 Sec.	Sound off
3	3	Cycle of RF channel check	2 Sec.	5 Sec.	10 Sec.	15 Sec.	30 Sec.	60 Sec.
4	4	Cycle of portable receivers	1 Sec.	2 Sec.	3 Sec.	5 Sec.	7 Sec.	10 Sec.
5	5	Relay on time after incoming Call	1 Sec.	2 Sec.	3 Sec.	5 Sec.	10 Sec.	15 Sec.
6	6	Frequency channel	1	2	3	4	-	-
7	1 + 2	Cycle of signalers	1 Sek.	2 Sek.	3 Sec.	5 Sec.	7 Sec.	10 Sec.
8	2 + 3	Signaler LED-Mode	Light	Blink	Flash			
9	3 + 4	Signaler alarm	Off	10 Sec.	30 Sec.	1 Min.	5 Min.	10 Min.
10	5 + 6	Configuration mode	-	-	-	-	-	-

Table 2: Parameter and selectable values (settings in bold are default values)

4.5.2 Volume Setting of Call Sound

The volume of the calling sound may be adapted to individual needs. Therefore open the battery housing (picture 1 - 3) take off the last battery and tune with a small screw driver until you have got the desired volume.

Picture 1: Back of Minicall-12

5 Additional functions

5.1 Battery supply

Pull the isolating stripe on the back of the housing (picture 1 - 2) in order to activate the battery supply which is only a back up supply if there is an AC-supply failure. The LED "DATA" stops flashing. Please take into account, that the used batteries only have a limited capacity when taking off the AC-adaptor. As long as the AC-adaptor is connected, the batteries will not be discharged. Battery supply is only to overcome an AC-Supply failure (i.e. changing the placement).

5.2 Relay Contact

The Minicall12 is equipped with a potential free relay contact which closes for a certain time (see table 2) each time a call from a transmitter is coming in.

Only devices with <40V are allowed to be switched with this contact.

You may switch on for example signalling lamps or phonetransmitters. Please take into account, that this feature is only a switching contact and the external device connected may require an external supply.

If you want to use this contact, you have to open the housing (picture 1 - 4) and use the connector indicated.

Picture 2: Tune volume of internal loudspeaker

6 Maintenance

The device is equipped with 5 alkaline batteries (AA) in order to overcome AC-supply failures for up to some hours (depending on state of the device). If batteries have to be changed LED „DATA“ flashes twice a second.

If this uninterruptible power supply is an important function for your application check regularly the battery state by taking off the AC-adaptor. The green LED "ON" flashes once a second if battery capacity is ok. If it flashes twice a second so you have to change the batteries.

Detailed instructions on installation and maintenance can be found in the corresponding documentation for the Minicall-12. We recommend a complete maintenance according to the SVS maintenance protocol at least once a year. System states and error messages have to be monitored constantly.

We recommend to change batteries every year. Therefore you have to open the battery housing (picture 1 -3). Only use alkaline 1.5 V size AA batteries and pay attention to the right polarity.

7 Technical Details

Antenna: integrated

Dimension: 120mm x 140mm x 35mm

Supply 1: AC adaptor 230V AC – 12V DC

RF frequency: 433.62 MHz (433.82 – 434.02 – 434.22 MHz optionally)

Output: isolated switch over contact (24V DC / 1 A)

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Declaration of Conformity in accordance with the Directive 2014/53/EU (RED)

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2014/53/EU European Radio Equipment Directive (RED)
2011/65/EU Restriction of Hazardous Substances (RoHS)

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SVS Nachrichtentechnik GmbH